

# MemoPart™

Patent Ductus Arteriosus

(PDA) Occluder CE

MemoPart™ PDA Occluder is a percutaneous, transcatheter device applied to the closure of patent ductus arteriosus (PDA). The retention skirt guarantees the secure positioning in the ampulla of the ductus.



Cylinder Shape  
PDA Occluder




Cone Shape  
PDA Occluder



Hubless Cylinder Shape  
PDA Occluder




Hubless Cone Shape  
PDA Occluder




MemoPart™ Cylinder Shape PDAO

Catalogue No.	Device Size	A Aortic Disc Diameter (mm)	B Aortic Waist Diameter (mm)	C Pulmonic Waist Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WBFDQ-I 04	04	8	4	4	4.0	6F
WBFDQ-I 05	05	9	5	5	5.0	6F
WBFDQ-I 06	06	10	6	6	6.0	6F
WBFDQ-I 07	07	11	7	7	6.5	6F
WBFDQ-I 08	08	12	8	8	6.5	6F
WBFDQ-I 09	09	13	9	9	7.0	6F
WBFDQ-I 10	10	14	10	10	7.5	6F
WBFDQ-I 11	11	15	11	11	8.0	7F
WBFDQ-I 12	12	16	12	12	8.5	7F
WBFDQ-I 13	13	17	13	13	8.5	8F
WBFDQ-I 14	14	18	14	14	9.5	8F
WBFDQ-I 16	16	21	16	16	10.5	9F
WBFDQ-I 18	18	23	18	18	10.5	10F
WBFDQ-I 20	20	25	20	20	12.0	12F
WBFDQ-I 22	22	27	22	22	12.0	12F




MemoPart™ Hubless Cylinder Shape PDAO

Catalogue No.	Device Size	A Aortic Disc Diameter (mm)	B Aortic Waist Diameter (mm)	C Pulmonic Waist Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WTWBFDQ-I 04	04	8	4	4	4.0	6F
WTWBFDQ-I 05	05	9	5	5	5.0	6F
WTWBFDQ-I 06	06	10	6	6	6.0	6F
WTWBFDQ-I 07	07	11	7	7	6.5	6F
WTWBFDQ-I 08	08	12	8	8	6.5	6F
WTWBFDQ-I 09	09	13	9	9	7.0	6F
WTWBFDQ-I 10	10	14	10	10	7.5	6F
WTWBFDQ-I 11	11	15	11	11	8.0	7F
WTWBFDQ-I 12	12	16	12	12	8.5	7F
WTWBFDQ-I 13	13	17	13	13	8.5	8F
WTWBFDQ-I 14	14	18	14	14	9.5	8F
WTWBFDQ-I 16	16	21	16	16	10.5	9F
WTWBFDQ-I 18	18	23	18	18	10.5	10F
WTWBFDQ-I 20	20	25	20	20	12.0	12F
WTWBFDQ-I 22	22	27	22	22	12.0	12F



MemoPart™ Cone Shape PDAO

Catalogue No.	Device Size	A Aortic Disc Diameter (mm)	B Aortic Waist Diameter (mm)	C Pulmonic Waist Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WBFdq-II 06	0406	10	6	4	6.0	6F
WBFdq-II 08	0608	12	8	6	6.5	6F
WBFdq-II 10	0810	14	10	8	7.5	6F
WBFdq-II 12	1012	16	12	10	8.5	7F
WBFdq-II 14	1214	18	14	12	9.5	8F
WBFdq-II 16	1416	20	16	14	10.5	9F
WBFdq-II 18	1618	23	18	16	10.5	10F
WBFdq-II 20	1820	25	20	18	12.0	12F
WBFdq-II 22	2022	27	22	20	12.0	12F



MemoPart™ Hubless Cone Shape PDAO

Catalogue No.	Device Size	A Aortic Disc Diameter (mm)	B Aortic Waist Diameter (mm)	C Pulmonic Waist Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WTWBFdq-II 06	0406	10	6	4	6.0	6F
WTWBFdq-II 08	0608	12	8	6	6.5	6F
WTWBFdq-II 10	0810	14	10	8	7.5	6F
WTWBFdq-II 12	1012	16	12	10	8.5	7F
WTWBFdq-II 14	1214	18	14	12	9.5	8F
WTWBFdq-II 16	1416	20	16	14	10.5	9F
WTWBFdq-II 18	1618	23	18	16	10.5	10F
WTWBFdq-II 20	1820	25	20	18	12.0	12F
WTWBFdq-II 22	2022	27	22	20	12.0	12F