

MemoPart™

Atrial Septal Defect

(ASD) Occluder CE

MemoPart™ ASD Occluder is a percutaneous, transcatheter device applied to the closure of atrial septal defects (ASD). It is also available for the patient who require a closure of fenestrated Fontan.



ASD Occluder




Hubless ASD Occluder



Multi-Fenestrated
ASD Occluder




Hubless Multi-Fenestrated
ASD Occluder



MemoPart™ ASDO


Catalogue No.	Device Size	A Waist Diameter (mm)	B LA Disc Diameter (mm)	C RA Disc Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
FQFDQ-1 06	06	6	16	14	3.0	8F
FQFDQ-1 07	07	7	21	17	3.0	8F
FQFDQ-1 08	08	8	18	16	3.0	8F
FQFDQ-1 09	09	9	23	19	3.0	8F
FQFDQ-1 10	10	10	20	18	3.0	8F
FQFDQ-1 11	11	11	25	21	3.0	9F
FQFDQ-1 12	12	12	22	20	3.0	9F
FQFDQ-1 13	13	13	27	23	3.0	9F
FQFDQ-1 14	14	14	24	22	3.0	9F
FQFDQ-1 15	15	15	29	25	4.0	9F
FQFDQ-1 16	16	16	30	26	4.0	9F
FQFDQ-1 17	17	17	31	27	4.0	9F
FQFDQ-1 18	18	18	32	28	4.0	9F
FQFDQ-1 19	19	19	33	29	4.0	10F
FQFDQ-1 20	20	20	34	30	4.0	10F
FQFDQ-1 22	22	22	36	32	4.0	10F
FQFDQ-1 24	24	24	38	34	4.0	10F
FQFDQ-1 26	26	26	40	36	4.0	12F
FQFDQ-1 28	28	28	42	38	4.0	12F
FQFDQ-1 30	30	30	44	40	4.0	12F
FQFDQ-1 32	32	32	47	42	4.0	12F
FQFDQ-1 34	34	34	49	44	4.0	14F
FQFDQ-1 36	36	36	51	46	4.0	14F
FQFDQ-1 38	38	38	54	50	4.5	14F
FQFDQ-1 40	40	40	56	52	4.5	14F
FQFDQ-1 42	42	42	58	54	4.5	14F
FQFDQ-1 44	44	44	60	56	4.5	14F
FQFDQ-1 46	46	46	62	58	4.5	14F
FQFDQ-1 48	48	48	64	60	4.5	14F
FQFDQ-1 50	50	50	66	62	4.5	14F



MemoPart™ Hubless ASDO


Catalogue No.	Device Size	A Waist Diameter (mm)	B LA Disc Diameter (mm)	C RA Disc Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WTFQFDQ-1 06	06	6	16	14	3.0	8F
WTFQFDQ-1 07	07	7	21	17	3.0	8F
WTFQFDQ-1 08	08	8	18	16	3.0	8F
WTFQFDQ-1 09	09	9	23	19	3.0	8F
WTFQFDQ-1 10	10	10	20	18	3.0	8F
WTFQFDQ-1 11	11	11	25	21	3.0	9F
WTFQFDQ-1 12	12	12	22	20	3.0	9F
WTFQFDQ-1 13	13	13	27	23	3.0	9F
WTFQFDQ-1 14	14	14	24	22	3.0	9F
WTFQFDQ-1 15	15	15	29	25	4.0	9F
WTFQFDQ-1 16	16	16	30	26	4.0	9F
WTFQFDQ-1 17	17	17	31	27	4.0	9F
WTFQFDQ-1 18	18	18	32	28	4.0	9F
WTFQFDQ-1 19	19	19	33	29	4.0	10F
WTFQFDQ-1 20	20	20	34	30	4.0	10F
WTFQFDQ-1 22	22	22	36	32	4.0	10F
WTFQFDQ-1 24	24	24	38	34	4.0	12F

Your **Needs** • Our **Design**



MemoPart™ Multi-Fenestrated ASDO

Catalogue No.	Device Size	A Waist Diameter (mm)	B LA Disc Diameter (mm)	C RA Disc Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
FQFDQ-II 06	06	6	30	22	3.0	9F
FQFDQ-II 08	08	8	32	24	3.0	9F
FQFDQ-II 10	10	10	34	26	3.0	10F
FQFDQ-II 12	12	12	36	28	3.0	10F
FQFDQ-II 14	14	14	38	30	3.0	10F
FQFDQ-II 16	16	16	40	32	4.0	12F
FQFDQ-II 18	18	18	42	34	4.0	12F
FQFDQ-II 20	20	20	44	36	4.0	12F
FQFDQ-II 22	22	22	46	38	4.0	12F
FQFDQ-II 24	24	24	48	40	4.0	14F



MemoPart™ Hubless Multi-Fenestrated ASDO

Catalogue No.	Device Size	A Waist Diameter (mm)	B LA Disc Diameter (mm)	C RA Disc Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WTFQFDQ-II 06	06	6	30	22	3.0	9F
WTFQFDQ-II 08	08	8	32	24	3.0	9F
WTFQFDQ-II 10	10	10	34	26	3.0	10F
WTFQFDQ-II 12	12	12	36	28	3.0	10F
WTFQFDQ-II 14	14	14	38	30	3.0	10F

MemoPart™ (ASD) Occluder