




MemoPart™

Patent Foramen Ovale

(PFO) Occluder CE


MemoPart™ PFO Occluder is a percutaneous, transcatheter device applied to the closure of patent foramen ovale (PFO) for patients with the history of stroke, or the transient ischemic attacks (TIAs).



PFO Occluder




Hubless PFO Occluder



MemoPart™ PFOO

Catalogue No.	Device Size	A Waist Diameter (mm)	B RA Disc Diameter (mm)	C LA Disc Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
LYKFDQ-I 1818	1818	3.5	18	18	6.0	9F
LYKFDQ-I 1824	1824	4.0	24	18	7.0	9F
LYKFDQ-I 2424	2424	4.0	24	24	7.0	9F
LYKFDQ-I 2228	2228	4.5	28	22	7.0	10F
LYKFDQ-I 2828	2828	4.5	28	28	7.0	10F
LYKFDQ-I 2534	2534	5.0	34	25	7.0	12F
LYKFDQ-I 3434	3434	5.0	34	34	7.0	12F



MemoPart™ Hubless PFOO

Catalogue No.	Device Size	A Waist Diameter (mm)	B RA Disc Diameter (mm)	C LA Disc Diameter (mm)	H Waist Length (mm)	Smallest Recommended Sheath Size
WTLYKFDQ-I 1818	1818	3.5	18	18	6.0	9F
WTLYKFDQ-I 1824	1824	4.0	24	18	7.0	9F
WTLYKFDQ-I 2424	2424	4.0	24	24	7.0	9F
WTLYKFDQ-I 2228	2228	4.5	28	22	7.0	10F
WTLYKFDQ-I 2828	2828	4.5	28	28	7.0	10F
WTLYKFDQ-I 2534	2534	5.0	34	25	7.0	12F
WTLYKFDQ-I 3434	3434	5.0	34	34	7.0	12F